

Products

SLEEPSTYLE (INSOMNIA TREATMENT)
MICRO-CURRENT STIMULATOR



SLEEPSTYLE (INSOMNIA TREATMENT)

MICRO-CURRENT STIMULATOR

The new **NS-SLEEPSTYLE-3** is a compact and portable device which can be connected to an Android or IOS mobile device to access the mobile application for treatment of insomnia allowing users to use the device wherever they go.

This device uses Cranial Electrotherapy Stimulation (CES), which is a form of a non-invasive brain stimulation that applies a small, pulsed electric current across a person's head to treat insomnia through electrodes placed on the ear lobes.



CLINICAL APPLICATIONS

- Insomnia
- Anxiety
- Depression

DOWNLOAD PDF

Technical Specifications

Hardware			
Size Dimensions	187 x 113 X 48 mm		
Power Supply	2 AA batteries (Dual-channel maximum output)		
Operation Duration	≥ 20 hours		
Current Intensity	0-600uA		
Current Setting Resolution	50uA		
Stimulation Frequency	0.5 Hz	1.5 Hz	100 Hz
Stimulation Time	10 seconds	10 minutes	20 seconds
	30 minutes	1 hour	Continuous
Load Impedance	100Ω – 15K		
Current Output Intensity Error	50uA – 100uA (± 10%)		

Application Scope

Clinical Application

- Insomnia
- Anxiety
- Depression


Application Department

- Neurology
- Psychiatry
- Family Health

Components

SLEEPSTYLE COMPONENTS

Main Unit



Ear Clip Cable and Felt



Conduction Solution



USB Charging Cable



Sleep is an indispensable physiological process of human beings, which helps the repair of brain cells, recovery and enhancement of memory. A good sleep can ease fatigue, recover energy and improve immunity which promotes growth and development.

SleepStyle applies micro-current electrical stimulation therapy to effectively stimulate the brain's cerebral cortex and regulate the body's functions which results to the natural balance status of the body:

- During the day, the body and mind is active or excited.
- During the night, the body and mind is inhibited.

With continuous treatment, the user's nervous system and body adapts and their sleep disorder can be relieved/eliminated.

System Key Features

- Non-drug treatment for Insomnia – does not cause dependency
- Improves Sleep quality
- Can be used for regular treatment
- Compact and Mobile hand-held device
- Can be connected to Android or IOS device
- Easy to set-up – hardware to mobile device
- Easy to use and operate – mobile application
- Automatic shut-down
- Adjustable current intensity to ensure user's comfort



Treatment Type	Mechanism of Action	Schedule	Side Effect
Drug Therapy	The drug/medication is absorbed by the human body which then promote drowsiness	Taken before going to bed	<ul style="list-style-type: none"> • Dependence • Tolerance/Resistance • Liver and Kidney problems
Sleep Style	Stimulates the cerebral cortex through micro-current to achieve natural balance. (Day = Active ; Night – Inhibited)	During the day (Not recommended to use within 4 hours before going to bed)	<ul style="list-style-type: none"> • No Side Effects

Got a query for Neuro Style Product or Services?

CONNECT WITH US

GET IN TOUCH